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Fax: +49 30 8104-1827

E-Mail: dan.hodoroaba@bam

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Test material and software for the performance check of X-ray spectrometers attached to a scanning electron microscope

EDS-TM001

BAM provides the test material EDS-TM001, whose spectrum is sensitive to disturbances in the operation of energy dispersive X-ray spectrometers (EDS). All the measurements necessary for the performance check of an EDS attached to a scanning electron microscope (SEM) in compliance with ISO 15632:2002 „Microbeam Analysis – Instrumental specification for energy dispersive X-ray spectrometers with semiconductor detectors“ can be carried out by using of this sample only: (i) determination of the spectrometer energy resolution (FWHM) in dependence on the energy, (ii) check of the calibration of the energy scale, and (iii) check of the contamination state of the detector (L/K ratio).

The test material consists of an approx. 10 µm thick hard coating containing the elements C, Al, Mn, Cu und Zr, deposited on a steel (100Cr6 or V2A) or silicon substrate (s. Fig. 1).

The evaluation of the measurements is in principle possible with the software provided by the spectrometer manufacturer. Nevertheless, the evaluation can be readily performed with the software “EDS detector test” offered optionally by BAM. The software requires a detector provided with a polymer thin window.


Fig. 1 EDS-TM001 (left) and 10 kV spectrum (right)

The FWHM's of C-K and Mn-K α calculated with the BAM programme are shown in the two upper diagrams in Fig. 2. They are to be compared with the specifications of the spectrometer manufacturer. The Mn-L α /Mn-K α ratio is an indicator of a possible detector icing.

The diagram on the left bottom of Fig. 2 shows the spectrometer resolution in dependence on the X-ray energy.

The diagram on the right bottom of Fig. 2 shows the status of the energy calibration.


Fig. 2 Result of the spectrometer check

Berlin, August 2008

Contact:

Dr. Vasile-Dan Hodoroaba

Division VI.4

phone: +49 30 8104-3144

fax: +49 30 8104-1827

email: dan.hodoroaba@bam